


ENRAF-NONIUS X-RAY DIFFRACTION



Pictured above: On the left is a triple sample Guinier powder camera, on the right a temperature programmable Guinier camera, shown mounted on a sealed tube X-ray generator.

FOR YOUR DIFFRACTION APPLICATIONS, OUR EQUIPMENT LINE INCLUDES:

- X-Ray Generators: Conventional Sealed Tube and Rotating Anode Types
- Ambient Temperature, Temperature Programmable Analytic Guinier Cameras
- CAD-4 Single Crystal Diffractometer, FAST Area Detector
- Specialty Cameras
- EXAFS Systems
- Much More



For Information Please Contact **ENRAF-NONIUS** Or Your Local **ENRAF-NONIUS** Agent

ENRAF-NONIUS
P.O. Box 483
2600 AL Delft, Holland
Telephone 015-569230
Telex 84438083

ENRAF-NONIUS
390 Central Avenue
Bohemia, NY 11716
Telephone (516) 589-2885
Telex 960250

Delft New York Solingen Paris London Copenhagen

PD27

CEA

REFLEX FILMS FOR X-RAY DIFFRACTION

REFLEX 15 — Single coated for highest resolution in focusing cameras
Size: 18 x 24cm, box of 25 sheets FW
35mm x 35m roll

REFLEX 25 — Double coated for highest speed in Debye Scherrer and
single crystal cameras
Size: 5" x 5", box of 25 sheets FW
5" x 7", box of 25 sheets FW
35mm x 35m roll

Available from:

Cea AB
Box 174
S-152 01 STRANGNAS
Sweden
Telex: 10937 CEAPHOT S

Cea America Corp.
545 Steamboat Road
GREENWICH, CT 06830
USA
Telex: 6819246 SID UW

Cea (UK) Limited
75 Kilburn Lane
LONDON W10 4AW
England
Telex: 291569 BLISH G

Cea (Deutschland) GmbH
Postfach 100343
D-4330 MULHFIM/RUHR 1
West Germany
Telex: 856678 CEA

PD28

Materials Testing and Analysis

Offering

X-Ray Diffraction
X-Ray Fluorescence
Optical Emission
Atomic Absorption
Mechanical Testing
Optical & Scanning Electron
Microscopy

(201) 575-8665

1275 Bloomfield Avenue, Fairfield, N.J. 07006
Call for brochures: prices available

Short Course X-RAY POWDER DIFFRACTOMETRY and SPECTROMETRY May 19 — 23, 1986

The University of Texas at Austin
Theory — Hands-on Lab Work

- Powder Diffraction
Qualitative — Quantitative
- Search — Match
Computer and Manual
- Indexing
- X-Ray Fluorescence Spectroscopy

Instructor: Professor Hugo Steinfink
(512) 471-5233


Limit 12: Early Enrollment encouraged
Fee: \$900.00

Inquiries or request brochure:
Continuing Engineering Studies
College of Engineering, ECJ 10.324
The University of Texas at Austin
Austin, TX 78712
Phone: 512-471-3506

PD29

PD30

PROOF OF PURCHASE



**Your IBM-PC will love
MICRO-ID, and so will you!**

MICRO-ID offers you
Search/Match software for
powder diffraction analysis.
JCPDS files available include
Minerals, Common Phases,
and more.

Contact us for more
information. Ask about our
special Educational Package
and other X-ray software.

M I D
MATERIALS DATA, INC.
P.O. Box 791
Livermore, CA 94550
(415) 449-1084

PD31

the

Can you quickly perform XRD on particles as small as 5-15 μm ?

1973
IR-100 AWARD

SMALL DIAMETER CAMERA INSERTS

28.6mm or 14.3mm
available for 57.3mm camera

Modified XRD System
(with camera)

\$4,050⁰⁰

Modification of your existing
Philips camera

\$545⁰⁰

For technical data sheets, contact: Ian M. Stewart

Walter C. McCrone Associates, Inc., specialists in ultramicro-analysis for over 30 years, routinely analyze small particles with drastically reduced exposure time and, in turn, significantly reduced background using our modified PHILIPS camera.

Inquiries concerning your specific requirements and rates are welcome!

Please call:
Dr. Robert Z. Muggli

McCrone group



2820 SOUTH MICHIGAN AVENUE • CHICAGO, ILLINOIS 60616 • 312/842-7100

PRECISION GLASSBLOWING *of Colorado*

ICP GLASS SPECIALISTS

WE

Fabricate and repair the ICP glass accessories of all major-source manufacturers.

- Torches-torch extensions
- Nebulizers-spray chambers
- Quartz shields-tips for Jarrell Ash torches

Have most torches in stock for immediate delivery.





Provide custom fabrication and prototype development services.

YOU

are guaranteed the finest in craftsmanship and service.


Can realize substantial savings - often up to 50%.

Are invited to call or write us for further information.



2901 S. SHOSHONE
ENGLEWOOD, COLORADO 80110
(303) 761-9410

PD32



THE MINERAL POWDER DIFFRACTION FILE GROUP DATA BOOK

- Now available in one bound volume, this 1983 publication is designed to aid research in mineralogy or mineral related studies by placing together those minerals that have SIMILAR DIFFRACTION PATTERNS.
- In total 1600 mineral species have been arranged in 130 groups of which 44 have been divided into subgroups.
- Users may gain entry into THE MINERAL POWDER DIFFRACTION FILE GROUP DATA BOOK via the alphabetical mineral species index found in the back.

for further
information write:

JCPDs

INTERNATIONAL CENTRE
FOR DIFFRACTION DATA
1601 Park Lane, Swarthmore, PA 19081
(215) 328-9400



PD33

PD34

Powder Diffraction, Vol. 1, No. 1, March 1986

COMPONENTS AND SYSTEMS FOR SEMICONDUCTOR AND MATERIALS CHARACTERIZATION...

BLAKE INDUSTRIES



Blake Industries is the leading manufacturer of x-ray components and systems that meet the requirements of modern materials and semiconductor technology.


We develop—and constantly evolve—instrumentation sensitive to new materials and to the demands imposed by smaller and smaller semiconductor circuits. Our product designs reflect today's major emphasis on detecting material defects, as well as the need for analyzing epitaxial films by measuring lattice mismatch between the epitaxy and the film. And our instrumentation is known as the

most flexible in handling a wide range of materials with all the necessary radiation safety.


Other leading companies, such as IBM, Bell Labs, Westinghouse, TRW and AT&T, are putting Blake Industries' technology to work for them. Isn't it time yours did?

BLAKE INDUSTRIES, INC.

660 Jerusalem Road
Scotch Plains, NJ 07076
(201) 233-7240
Telex: 7109979581



PD35



The 23rd Annual
X-Ray Clinic[®]
SUNYA

Register now for these popular one- and two-week intensive courses which teach everything from basics to the latest advances.

	<u>Fundamentals</u>	<u>Advances</u>
X-Ray Spectrometry	June 2-6	June 9-13, 1986
X-Ray Powder Diffraction	June 16-20	June 23-27, 1986

Tuition: 1 week = \$900.00; 2 weeks = \$1700.00
Payable in U.S. dollars drawn on a U.S. bank.
Checks payable to: X-Ray Clinic, SUNYA.

For further information and brochure:
Prof. Henry Chessin, X-Ray Clinic, SUNYA
Department of Physics
1400 Washington Avenue, Albany, NY 12222
(518) 442-4512 or 4513

PD36

Selected Powder
Diffraction Data for
FORENSIC MATERIALS

Now available in two volumes, this new publication provides the powder diffraction analyst with a completely up to date listing of over 3200 forensic patterns.

The DATA BOOK lists forensic data cards from Sets 1-32 of the Powder Diffraction File in card number sequence. Among other materials, drugs, explosives, poisons, minerals, metals and alloys are included.

The comprehensive SEARCH MANUAL contains a Hanawalt numerical section, a Fink numerical section, an alphabetical index and a formula index.

JCPDS

For further information write:
**INTERNATIONAL CENTRE
FOR DIFFRACTION DATA**
1601 Park Lane, Swarthmore, Pa. 19081
(215) 328-9400

PD37

POWDER DIFFRACTION FILE

STANDARD REFERENCE FOR POWDER DIFFRACTION ANALYSIS

▶ **SET 35** NOW
AVAILABLE

- DATA CARDS
- MICROFICHE
- MAGNETIC TAPES & DISKS
- MICRO-COMPUTER S/M
- ELECTRON DIFFRACTION S/M
- SEARCH MANUALS
- CRYSTAL DATA
- SPECIALIZED PUBLICATIONS

JCPDS®

For further information write:
**INTERNATIONAL CENTRE
FOR DIFFRACTION DATA**
1601 Park Lane, Swarthmore, Pa. 19081
U.S.A.
(215) 328-9400 • Cable: JCPDS

PD38

**LET POWDER DIFFRACTION
DELIVER A TARGET MARKET
TO YOU**

Our international readership spends in excess of \$100,000,000 for X-ray instrumentation, and contracts for more than \$30,000,000 worth of after sales and services each year.

Capture your share of this thriving market by reserving space in our upcoming issues.

Call 201-697-0343 or Write:

Gerri G. Hahn
Advertising Manager
Powder Diffraction
3 Deer Haven Lane
Newfoundland, NJ 07435
USA

April 1, 1986 is the closing date for the June 1986 issue

INTRODUCING

Rigaku's D/Max-B Series

... the most complete, self-contained x-ray diffraction system available and at a remarkably low price!

Rigaku's new D/Max-B Series utilizes a self-contained, high speed, high capacity micro computer providing all the capabilities of a modern x-ray diffractometer system:

- Simultaneous operation of two instruments
- Automatic detector parameter selection
- Computer aided alignment
- Self diagnostics
- Full color graphics
- Data processing programs including software for search/match, quantitative & qualitative analysis, texture, residual stress determination, particle size, lattice parameter & indexing, deconvolution, microdiffraction and small angle scattering.

Amazing... all of these features available in one single unit... and the price... remarkably affordable! For more descriptive literature or a demonstration, call or write Rigaku/USA today!



Rigaku/USA, Inc.
3 Electronics Ave.
Danvers, MA 01923 U.S.A.
Phone: (617) 777-2446



Explore the advantage of the world's most sophisticated x-ray diffraction system.

Philips' computerized APD, the world's most sophisticated x-ray diffraction system, provides optimum results, yet requires little more than the touch of a button.

The Philips APD is totally user oriented. By combining a flexible hardware configuration with the most comprehensive software library available,

it offers ease of operation and maximum system capability.

Interactive graphic operation guides the user step by step. In routine studies, a batch program automatically collects the data, then provides final analysis and a printout of identified phases or a calculation of concentrations. Unequaled accuracy in identification of phases is accomplished by a fast, automated search/match routine. A wide choice of algorithms pro-

vides for quantitative and profile analyses.

For more sophisticated applications, the APD software package includes programs for indexing powder patterns, isostructural searches aiding in solid solutions analysis, preferred orientation calculations, automated texture analysis, incorporation of user-generated patterns into a data base and calculation of lattice parameters to provide ultra-accurate 'd' spacings.

Over 80 years of x-ray experience on an international scale has resulted in the largest user base and a global network of sales, training and service personnel that meet our customers' needs worldwide. The Philips computerized APD x-ray diffraction system is a prime example of our unique ability to anticipate the technology of tomorrow. For an even greater discovery of the capabilities of the APD or information on an upgrade of your present system, talk with Philips today.

Contact: N.V. Philips Analytical X-ray, Lelyweg 1, 7602 EA Almelo, The Netherlands. Tel. (31) 5490-18291. (In U.S.A.) Philips Electronic Instruments, Inc., A North American Philips Company, Analytical X-ray Group, 85 McKee Drive, Mahwah, NJ 07430. Tel. (201) 529-3800.



**Analytical
Equipment**

PHILIPS