

Detection of Rat Macrophages, an Immunoperoxidase Method (DAB-NiCl₂ with Methyl Green Countestain)

Don C. Wilson, University of Washington Medical Center

I. Slide Preparation

- 1) Using microtome, cut 5 mm thick tissue sections.
- 2) Place tissue sections on SuperFrost/Plus Fisher brand slides.
- 3) Place slides in a 45° C drying oven overnight.
- 4) Deparaffinize slides in 3 changes of xyliene for 5 minutes each.
- 5) Rehydrate the slides in the following order:
 - a) 100% ethyl alcohol - 3 changes, 5 minutes total
 - b) 95% ethyl alcohol - 2 changes, 4 minutes total
 - c) 70% ethyl alcohol - 1 change, 1 minute
- 6) Rinse slides in distilled water.
- 7) Place slides in 3% hydrogen peroxide for 5 minutes (save 5 mL for DAB solution).
- 8) Rinse slides in distilled water.
- 9) Rinse slides in PBS twice. Microwave antigen retrieval:
- 10) Place slides in a plastic rack (no metal). Rack must always be full. Load with blank slides if necessary.
- 11) Place slides in 200 mL of 10 mM citrate buffer, cover with a lid loosely and place at edge of microwave tray.
- 12) Microwave at high power setting for 10 minutes.
NOTE: Do not let fluid level evaporate below level of tissue. Stop microwave every 3 minutes and add distilled water.
- 13) Remove slides from microwave and let sit in citrate at room temperature for a full 30 minutes.
- 14) Rinse slides in phosphate buffered saline (PBS) and one change of

PBS/Triton for a total of 5 minutes.

II. Antibody Application

- 15) Incubate specimens (humidified chamber, room temperature) with 10% bovine serum albumin (BSA) in PBS for 20 minutes.
Use sufficient reagent to cover the specimen (usually 50-75 mL).
- 16) Wash slides in 2 changes of PBS and 1 change of PBS with BSA/Triton for a total of 5 minutes.
- 17) Prepare primary antibodies by diluting:
 - a) Control antibody (Mouse IgG) to 0.4 µg/mL in 1% BSA + PBS,
 - b) (Serotec ED-1) anti-rat macrophage antibody 1:100 in 1% BSA + PBS.
- 18) Apply an adequate amount (usually about 75 mL) of the primary antibody to cover the entire tissue. Include at least 1 slide as negative control (apply control mouse IgG).
- 19) Incubate the slides for 60 minutes in a moist chamber at room temperature.
At this time, prepare the secondary antibodies, (biotinylated anti-mouse IgG diluted to 6 µg/mL), and streptavidin/peroxidase complex. These solutions should be made up in 1% BSA + PBS.
Thaw the appropriate number of DAB aliquots.
Warm 175 mL of 0.05 M Tris buffer in a beaker (foil covered) in a water bath at 37° C. Warm 200 mL of 0.05 M Tris buffer in a slide boat in a water bath at 37° C.
- 20) Wash slides in 1 change of PBS and 1 change of PBS with BSA/Triton for a total of 5 minutes.
- 21) Incubate specimens with 10% horse serum* in PBS for 20 minutes.
Use sufficient reagent to cover the entire specimen (50-75 mL).
- 22) Wash slides in 1 change PBS.
- 23) Apply the biotinylated secondary antibody (anti-mouse IgG) and incu-

Denton Vacuum Homepage - Netscape

File Edit View Go Window Help

Back Forward Reload Home Search Guide Print Security Stop Netscape

Go to: [HTTP://WWW.DENTONVACUUM.COM/](http://www.dentonvacuum.com/)

From Sunglasses to Satellites.

Corporate Info Denton News

Denton Products Key Contacts

Representatives Directions

If it involves Thin Film Coatings,
it involves Denton.

<http://www.dentonvacuum.com>

[Desk II Cold Sputtercoater](#) | [DV 502A Carbon Coater](#) | [DV 401 Desktop Carbon Coater](#)
[HI RES 100 High Resolution Chromium Coater](#) | [DCP-1 Critical Point Dryer](#)
[Desk II TSC Turbo Sputter Coater](#) | [BTT III Benchtop Carbon Coater](#)

DENTON VACUUM 1259 N. Church Street Moorestown, NJ 08057 Phone 609.439.9100 FAX 609.439.9111

bate for 30 minutes in a moist chamber at room temperature.

- 24) Rinse with 3 changes of PBS.
- 25) Cover sections with peroxidase labeled streptavidin (Avidin:Biotinylated enzyme Complex [ABC] solution) and incubate for 30 minutes in a moist chamber at room temperature.
- 26) Wash slides in 1 change of PBS and 1 change of 0.05 M Tris buffer (at 37° C in a water bath) for a total of 5 minutes.

III. Color Development


NOTE: DAB may be carcinogenic. Use rubber gloves and prepare under a fume hood. DAB is light sensitive and should be prepared in the dark and just prior to use.

- 27) We prepare a stock DAB (3,3'-diaminobenzidine (3,3',4,4'-tetraminobiphenyl) tetrahydrochloride) solution by mixing 5.0 g DAB in 132 mL Tris buffer. We then divide it into 4 mL aliquots, freeze them at -70° C, and use them to prepare the working DAB solution as described. The 4 mL aliquots of stock solution are enough to make 180 mL of the working solution which will stain 20 slides at a time.
- 28) Preparation of working DAB solution:
NOTE: Always add reagents in the following order:
 - a) Add 4 mL of thawed DAB to the warmed Tris buffer.
 - b) While stirring, add 1 mL of 8% NiCl₂.
 - c) Add 12 drops of 3% hydrogen peroxide with a Pasteur pipette.
 - d) Mix solution thoroughly.
- 29) Incubate slides in the DAB solution for 10 minutes.
- 30) Wash slides in distilled water for 1 change.
- 31) Counterstain slides in methyl green for 5 minutes.
- 32) Dehydrate slides in the following order:
 - a) 95% ethyl alcohol - 1 quick dip (methyl green is extremely soluble in water).
 - b) 100% ethyl alcohol - 2 changes for a few quick dips.
 - c) Leave the slides in the last change of 100% ethyl alcohol for 5 minutes.
- 33) Rinse slides in 3 changes of xylene for a total of 5 minutes.
- 34) Apply coverslips with Permount.


*Obtained from normal animals, heat treated, centrifuged and filtered through a .45 micron filter.

EDS Upgrades

Upgrade your existing EDX system with the world leader in affordable EDX upgrades.



Our new digital pulse processor provides increased throughput and will interface to most EDS detectors.



Precise Standardless quantitation for even complex matrices



"Active" digital imaging, linescans, and x-ray mapping.

Combining digital imaging with our EDS system produces a very simple yet powerful Microanalysis system.

IXRF has established it-self as a "leader" in the Microanalysis industry as most of the major EDS companies convert from analog to digital systems.

- Hundreds of systems installed
- Easy to use Windows 95 or NT
- Each system comes complete with Pentium PC
- Installation and training at your site
- New detectors available
- Morphology and Feature analysis
- Automatic Point Analysis
- Quantitative Mapping and Linescans

IXRF Systems, Inc.
15715 Brookford Dr.
Houston, TX 77059
Tel:281/286-6485 Fax:281/286/2660
Website: ixrfsystems.com

RMC

3450 South Broadmont Drive, Suite 100
 Tucson, AZ 85713
 (520)903-9366 Fax: (520)903-0132

Configure your next Ultramicrotome to fit your needs

. . . . The MTX

See all options at our web-site <http://www.rmc-scientific.com/microtomes/>
 Also visit us at the following trade shows and conventions:

ASCP Show: Oct. 18-21 1998, Washington, D.C.
Neuroscience Show: Nov. 8-11 1998, Los Angeles, CA
Veterinary Pathologists: Nov. 17-19 1998, St. Louis, MO
MRS Show: Nov. 30 - Dec. 4 1998, Boston, MA
ASCB Show: Dec. 12-16 1998, San Francisco, CA