

# Powder Diffraction

H. Hermann and M. Ermrich	Microabsorption Correction of X-Ray Intensities Diffracted by Multiphase Powder Specimens	189
L.S. Zevin and Sh.L. Zevin	Standardless Quantitative X-Ray Phase Analysis – Estimation of Precision	196
B.L. Davis, D.K. Smith and M.A. Holomany	Tables of Experimental Reference Intensity Ratios. Table No. 2, December 1989	201
T.C. Huang and P.W. Wang	X-Ray Powder Diffraction Analysis of a Photoconductor Material p-Diethylaminobenzaldehyde Diphenylamine Hydrazone	206
J.A. Emig and D.K. Smith	The Detection and Quantification of Low Concentrations of Quartz in Dolostone by X-Ray Powder Diffraction	209
E. Goldish	X-Ray Diffraction Analysis of Barium-Strontium Sulfate (Barite-Celestite) Solid Solutions	214
A. Borgna, B. Moraweck and P. Fessler	X-Ray Diffraction and EXAFS Study of Palladium Tetrammine Dichromate, $[\text{Pd}(\text{NH}_3)_4]\text{Cr}_2\text{O}_7$	217
F.N. Blanchard	X-Ray Powder Data for $\text{CaWO}_4$ , Synthetic Scheelite	220
Fang Zhou <i>et al.</i>	$(\text{Na}_{0.5}\text{Bi}_{0.5})\text{TiO}_3$ : Synthesis and X-Ray Powder Diffraction Data	223
J.M. Amigó <i>et al.</i>	Crystal Data for Nickel Cimetidine Chloride: $\text{Ni}(\text{CM})_2\text{Cl}_2 \cdot 2\text{H}_2\text{O}$	225
F.N. Blanchard	New X-Ray Powder Data for Gorceixite, $\text{BaAl}_3(\text{PO}_4)_2(\text{OH})_5 \cdot \text{H}_2\text{O}$ , an Evaluation of d-Spacings and Intensities, Pseudosymmetry and its influence on the Figure of Merit	227
D.F. Rendle <i>et al.</i>	X-Ray Powder Diffraction Data for 2-nitro-1-phenyl-prop-1-ene	231
J.L. Garin and J.A. Costamagna	X-Ray Powder Diffraction Data for Two Ammonium Copper Molybdates	233
	International Report	236
	Calendar of Meetings	237
	Meeting Reports	238
	Correspondents Reports	240
	Short Courses and Workshops	240
	Book Reviews	241
	Computer Comments	242
	Commercial Announcements	245
	Volume 4, 1989	Contents 247
		Author Index 249



Volume 4 Number 4 December 1989

Powder Diffraction An international journal of materials characterization



# SIEMENS

## Ask for performance - ask for the new D 5000 x-ray diffractometer

### Detector configurations include:

- PSD (Position Sensitive Detector).
- Solid-state detector for improved detector efficiency.
- Standard scintillation counter.

### Versatile goniometer design:

- Operates in horizontal or vertical mode for reflection or transmission measurements.
- Converts from  $\theta:2\theta$  geometry to  $\theta:\theta$  in your lab.
- Converts to parallel beam.


### Fully automated design:

- Microprocessor control of data collection.
- Host computer for experimental setup and data analysis.

### Complete set of attachments for specific applications:

- Eulerian cradle for texture.
- Open Eulerian cradle for texture/stress.
- Automatic sample changer.
- High and low temperature stages.
- Programmable sample rotation.

Available with computer-controlled variable aperture slits in incident and diffracted beam paths to control divergence and reduce background.



For maximum performance in your lab, Siemens has engineered a new generation of x-ray diffractometer designed to meet your requirements. The D 5000 features a new goniometer design manufactured to the strictest tolerances. Complete with application-specific attachments and advanced software routines, the D 5000 has the precision and flexibility to outperform other systems in every x-ray application.



### *Your Solution is Siemens*

In USA & Canada: contact Siemens Analytical X-Ray Instruments, Inc. • 6300 Enterprise Lane • Madison, WI 53719-1173 • (608) 276-3000 • FAX (608) 276-3015  
Worldwide Contact: Siemens AG, Analytical Systems E 689 • D 7500 Karlsruhe 21 • P.O. Box 21 1262 • Federal Republic of Germany • Tel. (0721) 595-4295



# SIEMENS

## The art of texture analysis



Siemens new 3D software is the elegant solution to preferred orientation problems in materials science.

**Siemens . . . the company you can believe in.**

Worldwide Contact: Siemens AG, Analytical Systems E 689 • D 7500 Karlsruhe 21 • P.O. Box 21 1262 • Federal Republic of Germany • Tel. (0721) 595-4295  
In USA & Canada contact Siemens Analytical X-Ray Instruments, Inc., a joint venture of Siemens & Nicolet.  
6300 Enterprise Lane • Madison, WI 53719-1173 • Tel (608) 276-3000 • FAX (608) 276-3015

PD-2



# SPEX Handbook Latest Edition Available Now . . .

## Handbook of Sample Preparation and Handling

The latest Edition of the SPEX Handbook answers questions about

- producing homogeneous powders
- making pressed or fused pellets
- controlling contamination effects
- eliminating particle size effects
- handling liquid samples
- selecting and ordering SPEX products

The SPEX Handbook is an indispensable resource for spectroscopists. And it's **FREE**. To get your copy, clip and mail the coupon below.

Or call us today at

**800-LAB-SPEX**  
**(800-522-7739)**

PD-3

Rush me a free copy of your handbook!

Name/Title \_\_\_\_\_

Institution/Company \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_

State \_\_\_\_\_ Zip \_\_\_\_\_

Phone \_\_\_\_\_ Application \_\_\_\_\_

**SPEX**

INDUSTRIES, INC.  
3880 PARK AVE. • EDISON, N.J. 08820  
TEL: 201-549-7144 • FAX: 201-603-9647

<https://doi.org/10.1017/S0885715600013683> Published online by Cambridge University Press

Powder Diffraction, Vol. 4, No. 4, December 1989



# Powder Diffraction

An International Journal of Materials Characterization

## Deane K. Smith

Editor in Chief

Dept. of Geosciences  
The Pennsylvania State University  
239 Deike Building  
University Park  
Pennsylvania 16802 USA

## Jan W. Visser

European Editor

Henry Dumontlaan 81  
2614 GL, Delft  
Netherlands

## Richard N. Rose

Assistant Editor and  
Manager of Publication

## Ron Jenkins

Managing Editor

JCPDS-International Centre  
for Diffraction Data  
1601 Park Lane  
Swarthmore  
Pennsylvania 19081 USA

## Brian H. O'Connor

Editor for Australia  
and New Zealand

Curtin University  
Dept. of Applied Physics  
GPO Box U 1987  
Perth 6001 Western Australia  
Australia

## Mary M. Rossi

Assistant to the  
Managing Editor

## Gregory J. McCarthy

Editor for New Diffraction Data

North Dakota State University  
Department of Chemistry  
Fargo, North Dakota 58105-5516 USA

## Helein D. Hitchcock

Departments Editor

NASA DM-MSL-1  
Kennedy Space Center  
Florida 32899, USA

## Eloise Humez Evans

Editorial Assistant

## Editorial Advisory Board

A. Albinati, Milan, Italy  
L.V. Azaroff, Storrs, Connecticut  
C. Baerlocher, Zurich, Switzerland  
C.S. Barrett, Denver, Colorado  
P. Bayliss, Calgary, Alta., Canada  
C.Z. Bojarski, Katowice, Poland  
A. Brown, Nykoping, Sweden  
L.D. Calvert, Melbourne, Australia  
D. Cox, Upton, New York  
W. Eysel, Heidelberg, West Germany  
J. Fiala, Plzen, Czechoslovakia  
V.A. Frank-Kamenetsky, Leningrad, U.S.S.R.  
L. Frevel, Midland, Michigan  
P. Gado, Budapest, Hungary  
H. Goebel, Munchen, West Germany  
Th. Hahn, Aachen, West Germany  
G.G. Johnson, Jr., State College, Pennsylvania  
Q. Johnson, Livermore, California  
J.I. Langford, Birmingham, U.K.

D. Louër, Rennes, France  
H.F. McMurdie, Washington, District of Columbia  
M.E. Mrose, Washington, District of Columbia  
M.H. Mueller, Argonne, Illinois  
M. Nichols, Livermore, California  
B.H. O'Connor, Bentley, Australia  
W. Parrish, San Jose, California  
B. Post, Brooklyn, New York  
E. Prince, I.U.Cr. Representative  
R. Shirley, Guildford, U.K.  
R. L. Snyder, Alfred, New York  
W. Su, Changchun, China  
P. Suortti, Helsinki, Finland  
J.W. Visser, Delft, Netherlands  
S. Weissmann, Piscataway, New Jersey  
A.J.C. Wilson, Cambridge, U.K.  
T. Yamanaka, Tokyo, Japan  
R.A. Young, Atlanta, Georgia  
L. Zevin, Beer-Sheva, Israel

## Publisher

JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, Pennsylvania 19081, U.S.A.

*Powder Diffraction* is a journal of practical technique, publishing articles relating to the widest range of application – from mineral analysis to epitactic growth of thin films and to the latest advances in software. Although practice will be emphasized, theory will not be neglected, especially as its discussion will relate to better understanding of technique.

*Powder Diffraction* is published four times annually by the JCPDS-International Centre for Diffraction Data.

*Manuscript submissions.* The Editors will consider all manuscripts received, but assume no responsibility regarding them. Materials will be returned only when accompanied by appropriate postage.

*Subscriptions.* The annual subscription rate in the United States and Canada is \$39.50; Library; \$75.00; other than U.S.A., Canada and the Far East, the annual subscription is \$55.00. Subscriptions to the Far East, including Japan, China, Taiwan, Malaysia, the Philippines, Indonesia and Korea should be made via Sanyo Information System Corp., Koho Bldg. 11-8, 1-Chome, Kayabacho Nihonbashi Chuo-ku, Tokyo 103, Japan. Airmail delivery available for subscribers outside U.S./Canada for an additional cost of \$35.00 per volume (4 issues).

Payment may be made in U.S. dollars by company check, MasterCard, Visa or international money order. Please address communications to the publisher's office.

*Advertising.* For advertising rates and schedules contact the Publisher's Office JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, PA 19081, Telephone (215) 328-9405.

*Reprints and permissions.* Contact the Publisher's office.


*Postal Information.* Powder Diffraction (ISSN 0885-7156) is published quarterly for \$39.50 a year (U.S. and Canada) by the JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, Pennsylvania 19081. JCPDS principal office: 1601 Park Lane, Swarthmore, Pa. 19081. Julian Messick, Jr., General Manager. © 1988 JCPDS-International Centre for Diffraction Data. Postmaster: Send address changes to JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, Pennsylvania 19081.

International CODEN Service (Intercode): PODIEZ

ISSN 0885-7156


Telex 847170

FAX 215 328 2503







THE CHOICE OF THE FUTURE  
**SIMULTANEOUS XRD**  
**OVER 120° - 2 THETA**  
**- FAST - ACCURATE - RELIABLE -**




CPS 120 MOUNTED ON TGS (TRANSMISSION AND REFLECTION MODE)



SILICON  $\lambda = \text{CuK}\alpha$



QUARTZ  $\text{SiO}_2$   
 MONOCHROMATOR : QUARTZ



FAST DIFFRACTOMETER

**ACCESSORIES :**

CRYSTAT - FURNACE - SAMPLE SPINNER -  
 AUTOMATIC SAMPLE CHANGER - and  
 X-RAY GENERATOR etc...

**APPLICATIONS :**

POWDER, STRESS, and TEXTURE WITH  
 MATCHING SOFTWARE

**inel**

12, avenue de Scandinavie  
 91953 LES ULIS CEDEX FRANCE  
 PHONE 33 (1) 69.86.13.30 - FAX 33 (1) 69.86 14 19  
 TELEX 603965

**ASK FOR MORE INFORMATION AND ADDRESS  
 OF OUR REPRESENTATIVE IN YOUR COUNTRY**

PD-4



# Powder diffraction is our specialty.

In fact, it's all we do.



That's why we can deliver all those technological advances others have been promising. And we give you the largest selection of software there is. That's also why Scintag users enjoy the kind of after-sale support that helps them realize substantial returns on their investment, year after year.

Take our solid-state detector as an example. It's routinely available on the PAD V system, and it works for both horizontal and vertical configurations. We were the first to offer time averaged and signal processing data collection schemes for routine analysis. No one else comes close to our true microstep data collection technology ( $1^\circ=3200$  steps). Our continuously variable radius goniometer remains unique.

Peruse our software offerings, and you'll see how we've translated our powder diffraction expertise into practical solutions. In packages that solve everyday problems and address specialized needs for metallographic, crystallographic and thin film applications.

To help our users apply this advanced technology to their analytical problems, Scintag maintains the best customer support program in the industry. Through a combination of applications assistance and fast service turnaround, our users are able to fully exploit the capabilities of our systems.

Contact the powder diffraction specialists today. Better yet, ask for a demonstration, and see for yourself.



**Scintag, Inc.**

3350 Scott Boulevard, Building 52  
Santa Clara, CA 95054  
(408) 748-8544 TLX 171132



# Finally. An affordable, easy-to-operate alternative to high cost diffractometer automation.

## The Databox.

### The acquisition package:

The Databox is a stepping motor driver and data collection system specifically designed to control your diffractometer, all in the space of a two-wide NIM module. Just talk to it from a computer terminal or your PC using an incredibly friendly command language, and all your data acquisition needs are satisfied.

### The analysis package:


You can also buy the Databox bundled with MDI's Peak Identification and Micro-ID Search/Match software to run on your IBM PC (or compatible), giving you a complete x-ray control, acquisition, and analysis system.

### The bottom line:

For well under \$5000, the Databox will fully automate your x-ray equipment. Add the analysis software, and the total cost is just over \$10,000.

Finally. A system with proven reliability and performance at a reasonable price. To learn more, contact us at:


Radix Instruments  
1019 Stratford Ave.  
South Pasadena, CA 91030  
(818) 441-5351





# Supper

## Complete Debye-Scherrer Powder Camera Systems Accessories, and Film




114.6mm Powder Camera with optional viewing stand, illuminator, and Gandolfi "Single Crystal Randomizer" which is interchangeable with standard powder sample holders

Two Sizes: 114.6mm camera, 57.3 camera  
Easy to load and unload in the darkroom, the cameras are precision manufactured with no screw-on components. The collimator and beamtrap are accurately secured in position magnetically and a removable sample holder permits insertion of the specimen sample outside the body of the camera.

Guinier Camera for transmission and reflection photographs is also available.

**Simplified Specimen Line-up and Darkroom Handling Procedures Offer the Ultimate in Convenience**

- Gandolfi "Single Crystal Randomizer"
- Sturdy Viewing Stand and Illuminator
- Precision Film Punch and Trimmer
- Universal Track and Tripod Mount
- Large Inventory of Capillary Tubes
- Large Inventory of Kodak 35mm Film
- Film Measuring Device and Illuminator



Reads accurately to 0.01mm and eliminates eyestrain

**Manufacturers of Fine X-Ray Diffraction Equipment For Over 45 Years**

*Charles Supper Company*  
INCORPORATED


15 TECH CIRCLE, NATICK, MA 01760, U.S.A.  
TELEPHONE: (508) 655-4610  
1-800-323-9645

PD7

# Search Match. And Beyond.



## 1. Pattern Analysis

- ◆ Peak finding
- ◆ Profile fitting
- ◆ Size/Strain analysis
- ◆ Fast, interactive, complete




## 2. Search/Match

- ◆ CD ROM
- ◆ Subfiles
- ◆ User files



## 3. Tools

- ◆ Pattern simulation
- ◆ Indexing
- ◆ Cell refinement
- ◆ Laue patterns
- ◆ Others



**Complete XRD Analysis Software For Your PC. MDI.**

Materials Data, Inc.  
Post Office Box 791  
Livermore, California 94550  
415/449/1084  
FAX 415/373/1659



PD8



# SEIFERT

## XRD 3000

### X-RAY DIFFRACTION SYSTEMS



XRD 3000 systems are required in fundamental research and industry for powder diffraction phase, texture and stress analysis, orientation determination of single crystals or investigations of thin film, amorphous and liquid substances. The modular architecture of hardware and software allows total flexibility in building a high performance, versatile system for the varying applications.

XRD 3000 systems can be arranged with MZ VI goniometer, C 3000 microstep controller, ID 3000 computer-controlled generator, IBM PS/2 computer series with DOS or XENIX operating systems and SEIFERT XDAL 3000 software library.


Please contact us for detailed information. We have been specialized producers for diffraction equipment and industrial X-ray systems for NDT for more than 60 years.



**RICH. SEIFERT & CO. GmbH & Co KG**  
Bogenstr. 41 · D-2070 Ahrensburg · Phone (0 41 02) 76 - 0 · Telex 2 189 826 · Telefax (0 41 02) 7 61 89  
Röntgenwerk · West - Germany

PD9





## Second Announcement...

The first international congress on X-ray Analytical Methods for Materials Analysis will be held in Honolulu, August 12-16, 1991, at the Hilton Hawaiian Village.

## Congress Theme...

The major thrust of this meeting will be related to the practical aspects involved in X-ray methods for materials analysis. This will be in keeping with the tradition of the Australian X-Ray Analytical Association (AXAA), the Denver X-Ray Conference and the X-Ray Chemical Analysis Group of the Japan Society of Analytical Chemistry.

To be discussed will be the use of X-ray methods based on Powder Diffraction, Fluorescence, Surface Analysis, Absorptiometry, Column Electron Diffraction and Thin Film Characterization by X-ray Diffraction, and Trace Analysis and Thin Film Characterization by X-Ray Fluorescence.

A two day pre-congress workshop program will be held at the University of Hawaii at Hilo, on August 8 and 9.

### **PICXAM is organized by:**

The Australian X-ray Analytical Association  
The Denver X-ray Conference  
X-ray Chemical Analysis Group/  
Japan Society for Analytical Chemistry

### **PICXAM is cosponsored by:**

The International Centre for Diffraction Data  
The University of Denver  
The University of Hawaii


# PICXAM

## Pacific-International Congress on X-Ray Analytical Methods



# The smallest, lightest, 3kW x-ray generator you can buy.

The MAX3100 x-ray generator will power virtually all anode grounded x-ray tubes commonly used in wavelength dispersive fluorescence and x-ray diffraction applications. Numerous built-in safety features prevent system damage and maximize x-ray tube life.



WEIGHT: 72 lbs., 32.7 kg.

## SPECIFICATIONS:

Maximum power output:	3 kW
Maximum output voltage:	-60 kv DC
Maximum emission current:	50 mA
Output connector:	Federal 3 pin, L&S tied together
Voltage stability:	Better than .01% for $\pm 10\%$ line voltage variation.
Emission Current Stability:	Better than .01% for $\pm 10\%$ line voltage variation.
Input Voltage:	200-240 VAC 3-phase 50 or 60 hz 5-wire Delta or WYE system
Filament voltage:	0 to 12 VAC, 72 watts maximum.
Ripple:	Maximum—.1% RMS Typical—.05% RMS at full load and maximum current
Efficiency:	Better than 80%.
Cooling:	Forced air cooled.
Warranty:	Two years—return to our factory.



*American Instruments*, division of inel inc.

50 CAMPUS PLAZA DRIVE, EDISON, NJ 08837  
Analytical X-Ray Equipment

Phone: (201) 417-0070 Fax: (201) 417-0430


**CALL FOR OUR FREE BROCHURE WITH COMPLETE SPECIFICATIONS AND DETAILS**

PD-11




# PC-PDF


## The Powder Diffraction File on CD-ROM



**Search** on key fields within the data base.



**Search** on combinations of fields.




**Display** results within seconds.

**PC-PDF**, a low cost personal computer based system, makes the entire Powder Diffraction File available on a single CD-ROM disk ... and, through use of optimum packing and access algorithms, displays results within seconds.

- Use with an IBM-XT or equivalent system.
- Annually updated disks will be made from the full data base and will include all additions and corrections.
- Space saving: One disk contains data equivalent to a file cabinet full of Data Cards.

You can now enhance your present system for search, retrieval and display. For additional descriptive material and ordering information, contact:


JCPDS-International Centre for Diffraction Data  
1601 Park Lane, Swarthmore, PA 19081-2389 USA  
(215)328-9400 Telex: 847170



PD-12



# Introducing the world's fastest, high resolution powder diffractometry system ...



The new PDS120 is a complete, integrated system that provides superior results quickly and easily. It is the product of close cooperation between Enraf-Nonius, the world-leader in single-crystal X-ray diffractometry, and Inel, the manufacturer of a unique curved position sensitive detector.

## **Unique, patented anode offers superior results**

The anode wire of traditional linear detectors has been replaced by a patented blade anode. Its mechanical stability guarantees improved resolution and a better signal-to-noise ratio, allowing fast simultaneous data acquisition over a large angular range. The standard system comprises:

## **The advanced CPS 120 detector**

- ★ excellent resolution, count rate and efficiency
- ★ wide, 120° aperture
- ★ ideal for many applications: rapid phase analysis, thin film, kinetic and phase transition studies, etc.

## **The field-proven FR590 X-ray generator**

- ★ reliable, compact and versatile
- ★ advanced HF inverter technology
- ★ processor controlled
- ★ 3 kW, highly stabilized
- ★ comprehensive safety features.

## **A simplified goniometer**

- ★ horizontal or vertical version
- ★ expanded application possibilities
- ★ optional high- and low-temperature attachments
- ★ optional sample changer.

## **User-friendly software**

- ★ for data transfer, calibration, peak search and analysis
- ★ additional packages for a wide range of applications.

## the PDS 120, from Enraf-Nonius

### **Worldwide user support**

Enraf-Nonius is a world-leader in the field of X-ray diffraction. It can offer the user a complete range of instruments and accessories including; X-ray generators, diffraction cameras, low- and high-temperature attachments and control equipment, and complete structure determination diffractometer systems.

**Ask for more information today!**




**B.V. Enraf-Nonius Delft**  
P.O. Box 483  
2600 AL Delft  
THE NETHERLANDS  
Tel: (015) 698500  
Fax: (015) 619574  
Tlx: 38083

**Enraf-Nonius**  
390 Central Avenue, Bohemia  
New York 11716  
U.S.A.  
Tel: (516) 589 2885  
Fax: (516) 589 2068  
Tlx: 2960250

PD-13





## Look into $\mu$ PDSM . . . It's more than accuracy. It's performance.

Performance makes  $\mu$ PDSM the one search/match software system that approaches universality. It's the result of many years of research into putting all the complexities of qualitative XRD analysis in a single system for your personal computer—and now we've added graphics with a power that matches the unparalleled performance that established  $\mu$ PDSM's reputation.

The utility of function and clarity of content in  $\mu$ PDSM's graphics go beyond merely attractive presentation, to give you a valuable addition to the analytic power  $\mu$ PDSM puts at your finger tips. With options that include fully integrated CD-ROM PDF-2 retrieval, direct instrument control, data acquisition, and diffractogram analysis,  $\mu$ PDSM

gives you a flexibility of application to match the power of its performance.

The universal applicability and continuing enhancement of  $\mu$ PDSM reflect the evolving technology of a company dedicated to software innovation. And you realize the results of Fein-Marquart's commitment to progress through a continuing update policy that keeps you at the state of the art.

Look into  $\mu$ PDSM. It takes you beyond accuracy into the next generation of search/match performance.



Fein-Marquart Associates, Inc.  
7215 York Road · Baltimore, MD 21212  
(301) 821-5980


PD-14

XIII



# A new class of diffractometer with X-RED

X RED... an X Ray Electronic Diffractometer based on a Curved Position Sensitive detector ranging 120 degrees with a resolution of better than 0.03 dg. and its associated electronic, X-ray generator, water refrigerator unit, all included in an independent X-ray proof cabinet. Of course linked to a PC Olivetti Computer with its fine graphic performance, to display in real time a diffraction pattern. And more, automatic sample exchanger for capillaries, liquid, high, low temperature attachment. X-RED - CAN SAVE YOU A LOT OF TIME IN DATA COLLECTION SPECIALLY IN THE FIELD OF ROUTINE POWDER AND AMORPHOUS MATERIALS ANALYSIS.



PD-15

38066 RIVA DEL GARDA - Zona Industriale Baltera (ITALY) - Tel. 0464/553426 - Telex 400278 RIVTOUR

## Advance Your PDF Skills

### Attend a JCPDS-ICDD Short Course and Workshop

A series of 3 day courses is being offered by JCPDS-ICDD in the methods of identification and characterization of crystalline substances through the use of powder diffraction data.


The purpose of these courses is to build proficiency of the user in the interpretation of experimental powder data, especially in the application of information provided in the *Powder Diffraction File*.

Hands on work sessions. Use of both manual and computer/search methods. Descriptive workbook with illustrative examples.

**February 6-8, 1990**      **Denver, CO**  
**Embassy Suites Denver Airport (Hotel)**

**April 17-19, 1990**      **Florida**

For further information contact:  
Ms. Josephine Felizzi  
JCPDS-ICDD  
1601 Park Lane, Swarthmore, PA 19081  
USA (215) 328-9400 Telex: 847170



PD-16

JCPDS



# SIETRONICS

## SieRay 112

### Automation for X-ray Diffractometers

(For a lot less than you thought)



- Fit a stepping motor to a manual powder diffraction goniometer, with standardised kits available for most models. Alternatively an existing stepping motor can be utilized.
- Use your existing counting electronics, or purchase the SieRay 112 with built-in counting electronics (as pictured above), AND a new detector if necessary.
- Speed up your data reduction with automatic background stripping, peak finding and on-screen editing. There's even a manual peak finder for difficult cases.
- Experience the power of having the JCPDS® data base on-line, via Fein-Marquart Associates "uPDSM"® or using "SIROQUANT"® pattern synthesis/stripping quantitative software.
- Extend your automation by controlling auto sample loaders or driving multiple goniometers from the one interface. Versions of 112 software are also available for other XRD microprocessors.
- Consider the SieRay 113 automation for x-ray spectrometers — just as flexible and versatile as the SieRay 112 and just as cost-effective.
- Installation and commissioning available — or do it yourself in 1 day for a SieRay 112.


**SIETRONICS PTY LIMITED**

P.O. Box 84, HAWKER, A.C.T., 2614 AUSTRALIA

Phone: (062) 51 6611 : Telex AA62754 : Fax: (062) 51 6659

PD-17

## X-RAY DIFFRACTION SERVICES



### Structure Determination

- Unit Cells
- Structure Solution
- Absolute Configuration

### Material Characterization

- Phase Identification
- Percent Crystallinity
- Polymorphism
- Hydrated/Solvated Forms

- Method Development
- QC/QA Testing
- JCPDS Search
- Quantitative Studies

### Support Services

- Quantitative Elemental Analysis
- Mass Spectrometry
- Scanning Electron Microscopy/EDX
- Auger Electron Spectroscopy/SIMS

**ORS** ONEIDA RESEARCH SERVICES, INC.

One Halsey Road  
Whitesboro, NY 13492  
Phone (315) 736-3050

PD-18

## X-ray Diffraction Reference Standards and Zero-background Sample Plates Custom Designed and Built for any Application

Your first step to improved x-ray diffraction results should be to contact The Gem Dugout for quality diffraction alignment standards and zero-background plates. And the next step is successful x-ray diffraction results.

**The Gem Dugout**  
1652 Princeton Drive  
State College, PA 16803  
(814) 865-5782

PD-19



# SEIFERT

## X-RAY TUBES



The SEIFERT diffraction X-ray tubes SN 60 with normal focus and SF 60 with fine focus are featured by a high primary intensity at a homogeneous intensity distribution, spectral purity, and a long lifetime.



**SEIFERT X-RAY CORPORATION**  
2551 Industry Lane, P.O. Box 294, Fairview Village, Pennsylvania 19409  
TEL: 215-539-4700 · TWX: 510-660-0568 · FAX: 215-539-6031



## Editorial

### Short Course on Advanced Powder Diffraction Methods

The Mineralogical Society of America held a short course on advanced powder diffraction methods in St. Louis, Missouri, November 4-5. The organizers were David L. Bish and Jeffrey E. Post who also edited the course notes which have been published by MSA. Other participants included R.C. Reynolds, Jr., R. Jenkins, R.L. Snyder, D.K. Smith, S.A. Howard, L.W. Finger and R.B. Von Dreele. The course included topics of considerable interest in research and materials characterization whose applications reached well beyond the mineralogical community.

Following the introductory material on the principles of diffraction, there was considerable discussion on instrument alignment and experimental procedures with modern equipment. Special topics included sample preparation, quantitative analysis, diffraction by small and disordered crystals, powder diffraction software, profile refinement methods, Rietveld refinement, synchrotron and neutron powder diffraction. The presentations emphasized the modern applications and advances in techniques and equipment. Experiments with neutron diffraction and synchrotron sources were also presented. Many of the presentations emphasized the use of sources with high intensities and digitized data traces.

With the increasing availability and widespread use of computer-controlled instrumentation and the availability of more intense sources, powder diffraction has expanded its role in materials analysis very rapidly in recent years. One good example of significant improvement is the use of the entire digitized diffraction trace for analysis rather than just the derived  $d$ s and  $I$ s. Applications to quantitative analysis, crystal structure refinement, and crystal perfection studies were emphasized in the MSA presentations. Use of the full digitized trace as raw or processed data will become increasingly important in the future. These early applications will be enhanced as new procedures are developed.

In addition to using the full trace as a basis for analysis, there has been considerable emphasis on the fitting of selected profiles to the digitized data to achieve the most accurate position and intensity information and in the unravelling of complex clusters of diffraction peaks. Not only have there been additional analytical profiles evaluated, but the computational methods in the fitting procedures have been optimized for rapid calculations. Profile fitting is an easier procedure than deconvolution, and the results are proving to be as useful in the determination of crystallinity parameters on imperfect materials.

One section of the MSA program, prepared by this editor, was a review of the available computer programs for analyzing diffraction information. This list contains over 100 programs, including many not discussed elsewhere in the course, along with the names of individuals to contact for copies. The programs are classified according to application, and similar programs are compared with respect to the type of algorithm used and the information produced. They are not compared for efficiency or ease of use. PC, mini and mainframe programs were considered. This list will be valuable to any diffractionist.

Although the course was designed and presented by mineralogists for mineralogists, the material is of general interest to any materials scientist. The course notes were prepared to be used as text material for advanced diffraction instruction or as a reference source for techniques and articles on general and special topics. This publication was issued in November 1989 and is available from the Mineralogical Society of America, 1625 I Street N.W., Suite 414, Washington, D.C. 20006, U.S.A. The cost of the publication is \$20.00 U.S.

Deane K. Smith  
Editor-in-Chief