


Using topography, thermal conductivity mapping and micro thermal analysis identifies layers in a polymer film. Courtesy of Duncan Price, Loughborough.



Topography and electrostatic force images combine to show misalignment of contacts and implants in SRAM.



When is an SPM

more than a microscope?

(When it has PPT.™)

PPT is Proximal Probe Technology.™ And it's only available in Scanning Probe Microscopes from ThermoMicroscopes.

Proximal Probe Technology gives you a whole lot more than just an image. It gives you the full range of proximity measurements you need to get accurate, precise and meaningful data in the engineering terms you recognize. Whether you want to measure friction, stiffness, adhesion, thermal conductivity, electrical or magnetic forces, or a host of other material properties, there's a ThermoMicroscopes system with PPT that's just right for the job. And for process characterization at the nanoscale level, nothing compares to the power PPT gives you.




Topography and phase images show area and percent coverage of Teflon® and lubricant in support matrix. Courtesy of Steve Pratt, Kodak.

REGION GROUP	AREA	
	μm ²	%
Total area	67.69	100
Silicone lubricant	3.101	4.6
Teflon particles	7.007	10.4

PPT is available on the full line of ThermoMicroscopes' instruments. Explorer™ for maximum versatility. CP Research,™ the SPM choice of research scientists demanding flexibility. AutoProbe M5™ for practical analytical and industrial applications. And Aurora,™ the world's leading NSOM system with exclusive, patented tuning fork technology.

To learn more about PPT and proximity measurement tools, and to get a free copy of "The Practical Guide to Scanning Probe Microscopy," call, fax or visit our web site today.



AutoProbe M5



THERMOMICROSCOPES
 Spark Scientific Instruments TOPOMETRIX

1171 Borregas Avenue, Sunnyvale, CA 94089
 Telephone 408.747.1600 Fax 408.747.1601

www.thermomicro.com

THE DIFFERENCE YOU CAN MEASURE™

As Easy As It Looks



PEO at 60°C with evidence of melting clearly shown after 2 minutes.

- ▶ **Based on a Simple PC Interface.**
- ▶ **Intuitive Operation with a Mouse Alone.**
- ▶ **Capable of Imaging Under Either Ambient or High Vacuum and at Both High and Low Temperatures.**
- ▶ **Fully Vibrationally Isolated.**
- ▶ **Capable of Wet or Dry Imaging.**
- ▶ **Variable Gas Pressure/Temperature.**

The new JSPM-4200 is a multipurpose turnkey Scanning Probe Microscope (SPM) capable of imaging using a myriad of different modes and

JEOL USA, Inc., 11 Dearborn Road, Peabody, MA 01960
Tel: 978/535-5900 • Fax: 978/536-2205
Email: spm@jeol.com • WEB: <http://www.jeol.com>

NEW - Multi-Environmental

SPM

Scanning Probe Microscope

accessories. Imaging modes include STM, STS (Scanning Tunneling Spectroscopy), CITS (Current Imaging Tunneling Microscopy), AFM, LFM (Lateral Force Microscopy), Non-contact Mode AFM, Discrete Contact Mode AFM, Lateral Modulation LFM, MFM (Magnetic Force Microscopy), and Force Modulation Microscopy.

If you would like to learn more about JEOL SPM's please give us a call and arrange for your own hands-on demonstration so that we can prove to you that it really is "As Easy As It looks".

JEOL