

ONLY MONOJECT® SAFETY SYRINGES HAVE A TRUE TRANSPORT POSITION

OUR PATENTED TWO-STEP LOCKING MECHANISM

ASSURES SAFE TRANSPORT AND SAFE DISPOSAL.

JUST SLIDE THE SHIELD FORWARD UNTIL YOU FEEL

THE POSITIVE STOP AND HEAR THE AUDIBLE SNAP.

IT'S NOW SECURE FOR SAFE TRANSPORT.

FOR SAFE DISPOSAL, SLIDE THE SHIELD FORWARD

AGAIN, AND, AFTER THE AUDIBLE SNAP, TWIST TO

LOCK PERMANENTLY. A SAFETY SHIELD WITH A

TRUE TRANSPORT POSITION . . . ONE MORE REASON

WHY HEALTHCARE PROFESSIONALS HAVE MADE

MONOJECT NO. 1. NOW AVAILABLE IN A FULL

ARRAY OF SIZES.

