

Volume rendering of the head of a snake scanned with 8 $\mu$ m resolution with left top corner virtually removed.



## SKYSCAN 1272

Performance Beyond  
Researchers' Imaginations




free QR-reader app  
is available on the AppStore

- More than 200 Megapixel (14450 x 14450 pixels) in every virtual slice through objects,
- More than 2600 such slices can be reconstructed after a single scan,
- Due to phase-contrast enhancement, object details down to 0.4 $\mu$ m can be detected,
- 20-100kV X-ray source, 6-position filter changer, integrated micropositioning stage,
- World's fastest hierarchical 3D reconstruction (InstaRecon®) with 20x...100x speed-up,
- 2D/3D image analysis, task lists, user plug-ins, surface and volume rendering,
- 16-position sample changer with auto adjustment of magnification and scanning protocol,
- The results can be exported to iPhone / iPad / Androids for 3D rendering by supplied software.

bruker-microct.com  
microtomography

Innovation with Integrity



# GV10X *Patented, ultra-high-efficiency downstream plasma cleaner*

**“Absolutely the most advanced electron microscope cleaner on the market...”**

**Removes hydrocarbons in-situ from SEMs, FIBs, TEMs and other analytical tools — and pre-cleans samples ex-situ**

### 10X Faster Hydrocarbon Removal

Exclusive 5 mTorr plasma and >20W power generates greater O<sup>0</sup> density with longer mfp — giving GV10X significantly higher ashing efficiency than competitive cleaners.

### Much Gentler Cleaning

Gently oxidizes hydrocarbons into molecules that are pumped free of the chamber. No specimen damage from energetic species, sputtering or heat.


### Faster Start, Total Control

Universal software works on any EM — no system integration needed. Start using GV10X in two minutes. Easy-to-use GUI gives operator full control of all functions.

### Uninterrupted Specimen Analysis

Want to clean a specimen during analysis? It's easy. Takes just 5 minutes for GV10X to cycle On/Clean/Off — not an hour!

Total-Control Graphical User Interface	W	Wr	HH	MM	SS
Set	50	..	0	30	0
Act	00	00	00	00	28
Lost	00	00	00	00	00
0.0V	6.8E+3	Pressure in Torr, mBar, or Pascal			



**ibss Group**  
*The World's Most Advanced EM Cleaning*

**Get more information: admin@ibssgroup.com • 1 (415) 566-5774 • www.ibssgroup.com**

**CAMBRIDGE**


New to the **Advances in Microscopy and Microanalysis** book series!

## Scanning Electron Microscopy for the Life Sciences

Heide Schatten

*University of Missouri, Columbia*

US\$120.00; Hb: 978-0-521-19599-7; 312 pp



Recent developments in scanning electron microscopy (SEM) have resulted in a wealth of new applications for cell and molecular biology, as well as related biological disciplines. It is now possible to analyze macromolecular complexes within their three-dimensional cellular microenvironment in near native states at high resolution, and to identify specific molecules and their structural and molecular interactions. New approaches include cryo-SEM applications and environmental SEM (ESEM), staining techniques and processing applications combining embedding and resin-extraction for imaging with high resolution SEM, and advances in immuno-labeling. With chapters written by experts, this guide gives an overview of SEM and sample processing for SEM, and highlights several advances in cell and molecular biology that greatly benefited from using conventional, cryo, immuno, and high-resolution SEM.

### About the series

The Press currently publishes the Microscopy and Microanalysis (MAM) journal in conjunction with the MSA, which reaches 4,000 microscopists and is affiliated with 12 international microscopy societies. The series would be a natural development from this journal, and will take a broad view of the discipline, covering topics from instrumentation to imaging, methodology and analysis across physical science, materials science, biology and medicine. Books commissioned for the series will range from advanced undergraduate textbooks through to research and practitioner oriented monographs for researchers. The series aims to produce a coherent source of material, encouraging the communication and exchange of ideas across these divergent fields, ensuring that the series appeals to a broad community in the physical and life sciences.

### Forthcoming titles in this series:

**Microscopic Nanocharacterization of Materials**  
by Michael Isaacson

**Energy Filtered Electron Microscopy and Electron Spectroscopy**  
by Richard Leapman

**Dynamic Transmission Electron Microscopy**  
by Nigel Browning, Thomas LaGrange, Bryan Reed, Henning Stahlberg, Bradley Siwick



[www.cambridge.org/us](http://www.cambridge.org/us)  
800.872.7423



**CAMBRIDGE**  
UNIVERSITY PRESS

## Isn't it about time you had your own Digital Microscope?

Now you can with the portable, affordable uScopeMXII!

The uScopeMXII is a small digital desktop microscope you can use in your workplace or home office. It captures images from standard glass slides and sends them to your PC.

You can interactively browse slides with full control of focus, image processing, and location. You can also scan regions of

interest creating fully focused image sets.

The industry-standard USB interface makes it simple to plug in and start capturing images. It easily interfaces with your desktop or laptop PC and allows you to view and capture slide images in a wide variety of environments.

The uScopeMXII is manufactured in the United States.

### Features and Benefits

- **Overview and Objective Cameras**  
The uScopeMXII has an objective camera for scanning and an overview camera for navigating.
- **Automatic Focus**  
Images are automatically focused using configurable focus algorithms.
- **Portability**  
At a weight of about 5 lbs., the uScopeMXII is highly portable.
- **Easy to Use**  
The uScope Navigator software simplifies scanning and browsing.
- **Self-Contained**  
The uScopeMXII is self-contained and includes the electronics, cameras, stage, and optics in a device about the size of a large external disk drive.
- **Full Imaging Control**  
User filters provide complete control over image processing and correction.



### Microscopes International, LLC

555 Republic Drive, Suite 119  
Plano, TX 75074-5498

Phone +1-214-785-2058  
FAX +1-214-785-2138  
Email [sales@usscopes.com](mailto:sales@usscopes.com)  
[support@usscopes.com](mailto:support@usscopes.com)

**Call us today for a quote  
or to schedule a demonstration**

# Amptek FAST SDD<sup>®</sup> for EDS (SEM) Applications




Transmission for Low Z Elements	
Be	13%
B	19.7%
Li	29%
C	43.9%
N	59.2%
O	62%
F	69%
Ne	72.9%
Na	75.1%
Mg	77.3%
Al	80.3%
Si	81.8%

Amptek's FAST SDD<sup>®</sup> detector for EDS use with SEMs utilizes new technology "C2 Series" X-ray windows (Si<sub>3</sub>N<sub>4</sub>) and has an excellent low energy response. Its high intrinsic efficiency makes it ideal for EDS XRF. See why Amptek detectors are the #1 choice of OEMs worldwide.

Water tight detector window!



**OEM's #1 Choice**



[www.amptek.com](http://www.amptek.com)

**AMETEK**  
MATERIALS ANALYSIS DIVISION

# Expand your Knowledge of Microscopy with MSA Membership!

Whether your primary focus is in optical, electron or scanning probe microscopy, the biological or the physical sciences, MSA takes your knowledge to the next level!

**Members Receive:**

- A personal subscription to MSA's official journal, *Microscopy and Microanalysis*, and MSA's popular bi-monthly magazine, *Microscopy Today*.
- Peer Networking through the Society's Focused Interest Groups and Local Affiliated Societies.
- Plus discounts on books, journals and other educational materials.
- MSA Awards Programs, Scholarships, Speaker Opportunities, and much more!



**Join MSA Today!**

For more information: visit [www.microscopy.org](http://www.microscopy.org) or call 1-800-538-3672




Visit the website often for details on Symposia, Exhibits, Biological and Physical Sciences Tutorials, In-Week Intensive Workshops, Posters, Sunday Short Courses and Pre-Meeting Congress!

# MICROSCOPY & MICROANALYSIS

August 2-6 • Portland, OR



<http://microscopy.org/MandM/2015>




Program Information | Information for Students | Meeting Awards | Exhibitor List & Floor Plan  
Online Registration (opens February 2015) | Portland Hotel Information (reservations available February 2015)



# SIA

**1 TO 50 MEGAPIXELS** live and slow scan  
**MAGNIFICATION FACTOR OF 1** on bottom mounted cameras  
**DIFFRACTION BEAM STOP** on side mounted cameras



Affordable TEM camera systems for research, education, healthcare, and industry since 2001

Scientific Instruments and Applications  
2773 Heath Lane • Duluth, GA 30096  
(770) 232 7785 • [www.sia-cam.com](http://www.sia-cam.com)

CAMBRIDGE JOURNALS


## Top Journal in Microscopy

2012 Impact Factor: 2.495\*

Published for the Microscopy Society of America

EDITOR: **Robert L. Price**  
University of South Carolina, USA

The only journal owned by scientists and published for scientists, *Microscopy and Microanalysis* provides original research papers in the fields of microscopy, imaging, and compositional analysis. This distinguished international forum is intended for microscopists in both biology and materials science.



Online submission at

<http://mc.manuscriptcentral.com/mam>

<http://journals.cambridge.org/mam>



\*2012 Journals Citation Report® Science Edition (Thomson Reuters, 2012)

Expand your Knowledge of Microscopy with

## MSA Membership!

Whether your primary focus is in the biological or the physical sciences, MSA takes your knowledge to the next level!

**Members Receive:**

- A personal subscription to MSA's official journal, *Microscopy and Microanalysis*, and MSA's popular bi-monthly magazine, *Microscopy Today*.
- Peer Networking through the Society's Focused Interest Groups and Local Affiliated Societies.
- MSA Awards Programs, Scholarships, Speaker Opportunities, and much more!

**Join MSA Today!**

For more information:  
visit [www.microscopy.org](http://www.microscopy.org)



or call 1-800-538-3672

// PERFORMANCE  
MADE BY ZEISS

# Your FIB-SEM for High Throughput Nanotomography and Nanofabrication **ZEISS Crossbeam 340 and 540**



Highest 3D  
Resolution



## Enjoy High Productivity with an Open 3D Nano-Workstation

With ZEISS Crossbeam 340/540 you experience an outstanding combination of imaging performance and analytical power. Process and prepare any sample on a nanoscopic scale with the next-generation FIB. Observe your process in parallel with the superb optics of the GEMINI column. Enjoy demanding experiments with your Crossbeam especially designed for long-term stability. Find out more and book a hands-on demonstration in one of our ZEISS Microscopy Labs now.



[www.zeiss.com/crossbeam](http://www.zeiss.com/crossbeam)



We make it visible.

June 7 - 12, 2015

# LEHIGH MICROSCOPY SCHOOL



Lehigh University, Bethlehem, PA USA

## MAIN COURSES

SCANNING ELECTRON MICROSCOPY  
AND X-RAY MICROANALYSIS  
**June 8-12**

INTRODUCTION TO SEM AND EDS  
FOR THE NEW OPERATOR  
**June 7**

## SPECIALIZED COURSES

FOCUSED ION BEAM (FIB):  
Instrumentation and Applications  
**June 8-12**

PROBLEM SOLVING:  
Interpretation and Analysis  
of SEM/EDS/EBSD Data  
**June 8-12**

QUANTITATIVE X-RAY  
MICROANALYSIS:  
Problem Solving using EDS  
and WDS Techniques  
**June 8-12**


SCANNING TRANSMISSION  
ELECTRON MICROSCOPY:  
From Fundamentals to  
Advanced Applications  
**June 8-12**

**For more information, contact:**  
Sharon Coe | 610.758.5133 |  
sharon.coe@lehigh.edu


Register and pay in full by  
April 14 to receive an early  
bird discount.

**45 YEARS OF EXCELLENCE**







4 channel simultaneous HRSTEM imaging of SrTiO<sub>3</sub> using 4 STEM detectors.




Combined 3D EDS map: Ni, Al, Ti, Mo, Cr, and Co.




Temperature driven Al aggregation in solar cell. Sample courtesy of Dr. S. Kraschewski, U. Erlangen.



HRSTEM of Si (110) at 200kV.



3D EDS maps at different angles.



HRTEM image of SrTiO<sub>3</sub> with Ceta 16M camera.

# Fast 2D & 3D Imaging and Chemical Analysis with Talos TEM

Whether you work with metals, light elements, or nanoparticles—from macro to atomic scales—Talos gives you insight into structure and composition for a complete understanding of material properties.

**High speed 2D and 3D EDS analysis:** Fast, multichannel, automated EDS acquisition and precise quantification to obtain multidimensional compositional maps.

**Fast navigation and imaging:** Highest throughput, highest resolution uncorrected S/TEM images—even for beam sensitive and challenging materials.

**Space to do more:** FEI TEMs have the space to accommodate additional *in situ* holders for dynamic imaging, diffraction, and tomography applications.



Watch the movie at [FEI.com/DiscoverTalos](https://www.fei.com/DiscoverTalos)



Explore. Discover. Resolve.



# THE NEW DESK V

## THE GOLD STANDARD IN SAMPLE PREPARATION



Now upgraded with storage for up to ten recipes and a color touch screen display, the Desk V HP delivers breakthrough electron microscopy sample preparation.

**DENTON VACUUM**  
BARRIERS BECOME BREAKTHROUGHS

Visit us at:  
[www.dentonvacuum.com/mt](http://www.dentonvacuum.com/mt)




---

### INDEX TO ADVERTISERS

---

Amptek	A18
Bruker MicroCT	A15
Bruker Nano GmbH	A3
Carl Zeiss Microscopy, LLC	A5, A21
Denton Vacuum	A24
EDAX, Inc.	A9
EMS/Diatome	A10, A11, A12, A13
FEI Company	A6, A23
Gatan	Cover 3
IBSS Group	A16
JEOL	A14, Cover 4
Microscopes International	A17
MSA Ads	A16, A18, A19, A20, A22
SI&A Scientific Applications	A20
TESCAN	Cover 2

---

# *In-situ* microscopy – uncompromised

There is always more to see when studying high speed reactions. That's why at Gatan we build large format cameras with unparalleled temporal resolution to push the boundaries of *in-situ* electron microscopy.

## OneView™

Use only ONE camera to capture high quality 16 megapixel still images AND video in all of your TEM applications.

- Always have a “live” experience with 25 fps at full 4k x 4k resolution
- > 300 fps at 512 x 512 resolution



[www.gatan.com/oneview](http://www.gatan.com/oneview)



[www.gatan.com/K2IS](http://www.gatan.com/K2IS)

## K2® IS



Resolve dynamic details *in-situ* with the fastest 4k x 4k direct detection camera available.

- Frame rates up to 1600 fps – sub-ms time resolution
- Highest available DQE results in images with lowest noise – every image is useful



# TAKE A FRESH APPROACH TO SEM

*Intuitive multi-touch screen interface • High spatial resolution • Unprecedented analytical versatility  
Easy navigation and operation • Fast imaging and data acquisition*



*Bring a whole new SEM experience  
to your lab with our InTouchScope  
series – designed to fit your budget  
and research applications.*

**JEOL**  
Solutions for Innovation

www.jeolusa.com  
salesinfo@jeol.com • 978-535-5900

[www.jeolusa.com/fresh](http://www.jeolusa.com/fresh)

