

INDEX TO ADVERTISERS

Diatome	Cover 2
EDAX	A9
Electron Microscopy Sciences, EMS	Cover 3, A7, A10
Hitachi	Cover 4, A8
Physik Instrumente	A4, A12
Tescan	A11

Seeing the future with clarity



The Clarity™ analysis system uses next-generation imaging technology for the direct detection of EBSD patterns. This revolutionary approach provides unparalleled pattern quality and distortion-free imaging, opening new doors to the evolution of EBSD pattern analysis.

Ultimate sharpness for maximum detail

Direct electron detection for zero distortion

True quantitative intensity measurements

Ideal for extremely beam sensitive materials

EDAX
EDAX logo with a red bar chart graphic below the text.

edax.com

AMETEK
AMETEK logo with a red triangle graphic to the left of the text.

your source
for cryo-SEM/
cryo FIB/SEM
preparation
equipment...

PP3010 Cryo-SEM/Cryo-FIB/SEM Preparation System

A highly automated, easy to use, column-mounted, gas-cooled cryo preparation system suitable for most makes and models of SEM, FE-SEM and FIB/SEM..

PP3004 QuickLok

Ambient temperature airlock for SEM, FIB/SEM, beamline and vacuum platforms.

PP3005 SEMCool

Non-airlock cryo cooling for SEM, FIB/SEM, beamline and vacuum platforms

PP3006 CoolLok

Cryo transfer systems for SEM, FIB/SEM, beamline and vacuum platforms

EMS has it!




PP3010 Cryo-SEM/Cryo-FIB/SEM Preparation System



PP3004 QuickLok



PP3005 SEMCool



PP3006 CoolLok



**PLEASE CONTACT US
FOR MORE INFORMATION...**

OUR MAIN INTERACTIVE WEBSITE:

www.emsdiasum.com

**Electron
Microscopy
Sciences**


P.O. Box 550 • 1560 Industry Rd.
Hatfield, Pa 19440

Tel: (215) 412-8400 • Fax: (215) 412-8450

email: info@emsdiasum.com

or stacie@ems-secure.com

Take gallium-free microsample preparation further



Array of micro-compression test pillars in UFG aluminium alloy

TESCAN AMBER X

- ✓ High throughput, large area FIB milling up to 1 mm
- ✓ Ga-free microsample preparation
- ✓ Ultra-high resolution, field-free FE-SEM imaging and analysis
- ✓ In-column SE and BSE detection
- ✓ Spot optimization for high-throughput, multi-modal FIB-SEM tomography
- ✓ Superior field of view for easy navigation
- ✓ Essence™ easy-to-use, modular graphical user interface



For more information visit

www.tescan.com

TESCAN AMBER X is based on the S8000 platform.

Nano-Positioning for Microscopy & Imaging



Piezo stages & positioners are essential tools for high-resolution microscopy, such as Super Resolution Microscopy or AFM. Their sub-atomic resolution and extremely fast response allow researchers to create higher-quality images faster. PI provides a large variety of fast positioning stages and piezo objective nano-positioners for 3D imaging (Z-stack acquisition), fast-focusing applications, and light sheet microscopy.



PIFOC® nano-focusing drive



*6-axis closed loop piezo stage for AFM:
highly linear, < .1nm resolution*




*Affordable XY & XYZ Piezo Stages
for SR Microscopy: P-545 Plnano®*

Based on 40+ years of experience with electromagnetic and piezo-ceramic motors, PI can quickly provide a solution for your precision positioning and automation projects. Design centers are located in the USA, Europe, and Asia with global support of 1300+ employees in 13 countries.

PI

Physik Instrumente
www.pi-usa.us
508-832-3456
AskPI@pi-usa.us



PRECISION | SPEED | STABILITY - MOTION CONTROL & POSITIONING SOLUTIONS

your source for
cutting-edge sample
preparation equipment...

GloQube® Plus Glow Discharge System for TEM Grids and Surface Modification

A cost-effective, compact and easy-to-use glow discharge system, designed to modify the surface of TEM grids to meet the requirements for imaging a variety of macromolecules.

Q150V Plus Automatic Coater

Optimized for high-vacuum applications, with an ultimate vacuum of 1×10^{-6} mbar. Use with a wide-range Penning/Pirani gauge enables the sputtering of oxidizing metals with ultra-fine grain sizes.

Q300T T Plus Triple Target Large Chamber, Turbo-Pumped Sputter Coater


Ideally suited for sputtering a single large diameter specimen up to 8"/200 mm or multiple smaller specimens over a similar diameter. Ideal for thin-film applications and SEM/FE-SEM.

Q300T D Plus Dual Target Sequential Sputtering System

Suitable for multi-layer sequential sputtering of two materials, with two independent sputtering heads, which allow sequential sputtering of two metals without the need to break vacuum.



GloQube® Plus



EMS has it!



Q300T D Plus



Q300T T Plus



Q150V Plus



**PLEASE CONTACT US
FOR MORE INFORMATION...**

OUR MAIN INTERACTIVE WEBSITE:

www.emsdiasum.com


**Electron
Microscopy
Sciences**

P.O. Box 550 • 1560 Industry Rd.
Hatfield, Pa 19440
Tel: (215) 412-8400 • Fax: (215) 412-8450
email: info@emsdiasum.com
or stacie@ems-secure.com


The Hitachi ArBlade 5000 Advanced Ion-Milling System

The Hitachi ArBlade 5000 Ion-Milling System is a highly advanced broad ion-beam system. The ArBlade 5000 equipped with a fast-milling Ar ion gun will increase your throughput and put your lab on the cutting edge.

**PUT YOUR LAB
ON THE CUTTING EDGE**



- 1 *New-generation hybrid instrument with cross-section milling and Flatmilling™*
- 2 *≥ 1,000 μm per hour cross-section milling on Si—Cut your milling time in half!*
- 3 *Cross-section widths up to 8 mm—Mill out massive areas with the new cross-section holder*



For more information: microscopy@hitachi-hta.com

Innovation ♦ Synergy ♦ Solutions