

Atomic Resolution at 50 – 300 kV Obtained Using Low Dose Rate HRTEM

B. Barton, C. Song and C. Kisielowski

National Center for Electron Microscopy, Lawrence Berkeley Nat'l Laboratory, Berkeley, CA

State-of-the-art high resolution TEM (HRTEM) using aberration correction allows for imaging of suitable samples with sub-Å lateral resolution [1]. The object exit wave can be recovered from through focal series [2]. In addition, the use of efficient, ultra-fast direct electron detectors [3] can enable the tracking of sample dynamics, such as atomic diffusion paths [4] and catalytic reactions [5]. Obtaining sufficient image signal-to-noise ratios (SNRs) for sub-atomic resolution generally depends on system stability and high-intensity illumination conditions [1].

However, it is proposed that sub-Å structural resolution in HRTEM is not only limited by instrumental parameters, but also by beam-induced atom movements such as the excitation of 3D lattice phonons [6]. Preliminary experiments with STO strongly suggest that exposure dose rates $<10 \text{ e}^-/\text{Å}^2\text{s}$ will significantly reduce such movements and allow for a more precise evaluation of exit wave phases.

In addition, a growing number of “soft” materials of interest [7], composed predominantly of light atoms (H,C,N,O...) do not permit HRTEM analysis with atomic resolution at present. Those samples are often beam-sensitive and subject to fast degradation under usual HRTEM conditions with dose rates of $> 1,000 \text{ e}^-/\text{Å}^2\text{s}$ (HDR conditions). In order to obtain atomically resolved structural information from soft materials, it is desirable to find imaging conditions with low dose rates of $\sim 10 \text{ e}^-/\text{Å}^2\text{s}$ (LDR conditions). Furthermore, displacement damage can be reduced by lowering the electron energy towards 50-20 keV. Therefore, combining LDR and low voltage conditions appears to be a promising approach to HRTEM of soft and hard materials.

We show that using the aberration-corrected TEAM 0.5 instrument equipped with a coherent high brightness electron source (XFEG) and a monochromator permits imaging of Si crystals at atomic resolution using 300 to 50 keV. Si(110) dumbbells can be resolved in both conventional HDR and LDR conditions down to a dose rate of $22 \text{ e}^-/\text{Å}^2\text{s}$. Aligned averages of 5 consecutive LDR images show sufficient SNR for exit wave reconstruction using the Gerchberg-Saxton algorithm. Comparison of exit wave phases from HDR and LDR datasets point towards a sharper phase signal for dumbbell atoms under LDR conditions. In the future, low-voltage LDR HRTEM, combined with ultra-fast direct detection and a physical phase plate, will be used to image hard-soft matter interfaces in materials relevant to chemical catalysis and harvesting of solar energy [8].

References

- [1] Kisielowski et al., *Microscopy & Microanalysis* 14 (2008)
- [2] Allen et al., *Ultramicroscopy* 100 (2009)
- [3] Bataglia et al., *Nuclear Instruments & Methods In Physics Research Section A* 608 (2009)
- [4] Alloyeau et al., *Phys. Rev. B* 80 (2009)
- [5] Kisielowski et al., *Angewandte Chemie, International Edition* 49, 2708-2710, 2010
- [6] Specht et al., *ChemCatChem*. 2011, in review
- [7] Nam et al., *Nature Materials* 9 (2010)
- [8] This investigation is part of the Soft Matter program, which is supported by the U.S.

Department of Energy under Contract No. DE-AC02-05CH11231. Electron microscopy at the NCEM is supported by the Office of Science, Office of Basic Energy Sciences, Scientific User Facilities Division, of the U.S. Department of Energy under Contract No. DE-AC02—05CH11231.




FIG. 1. : (a) Model of a Si crystal viewed along $\langle 110 \rangle$ direction. The characteristic dumbbells have a projected distance of 136 pm. (b) TEM image of the edge of a Si[110] plan view specimen obtained from a wafer by dimpling and subsequent etching with hydrofluoric acid (HF) and potassium hydroxide (KOH). (c) The etching process produces a “rooftop” topography favorable for HRTEM imaging.




FIG. 2. Si(110) “rooftop” sample imaged using the aberration-corrected, monochromated TEAM 0.5 microscope under conventional HDR and LDR conditions at 50, 80 and 300 keV electron energy. While individual LDR images show a low SNR, sufficient signal can be obtained by aligned averaging of 5 LDR images. 136 pm dumbbells can be resolved in HDR and LDR between 50 and 300 keV electron energy, even with a dose rate as low as $22 \text{ e}^-/\text{\AA}^2\text{s}$.




FIG. 3. Phase of exit wave reconstruction from HDR (a; $11,000 \text{ e}^-/\text{\AA}^2\text{s}$) and LDR focal series (b; $57 \text{ e}^-/\text{\AA}^2\text{s}$). (c,d) LDR conditions may yield sharper atomic phase images and better dumbbell peak separation.