



# Cambridge Journals Digital Archive


- ▶ 1770 – 2011
- ▶ Over 220 titles
- ▶ From Volume 1



# You did *what* in 5 minutes?



**3 LOG<sup>10</sup>**  
*C.diff* Reduction  
Total Room



**0.5 LOG<sup>10</sup>**  
*C.diff* Reduction  
4ft Radius Only\*

# 6x

Tru-D is more effective than pulsed *xenon's* 5 minute claim.

Every high-touch surface in the room requires disinfection. Tru-D performs at 99.9% efficiency guaranteeing disinfection of all surfaces including shadowed areas within any room from a single placement.

\* "Evaluation of a Pulsed Xenon Ultraviolet Disinfection System for Reduction of Healthcare-Associate Pathogens in Hospital Rooms," Infection Control & Hospital Epidemiology, Nerandzic, et. al, 2015