


**VACUTAINER®**  
QuikHeel™  
Safety Lancet



**VACUTAINER®**  
Safety-Lok™  
Needle Holder




**B-D®** Saf-T E-Z Set®  
Winged Infusion Set



**B-D®** Safety-Lok™  
Syringes



**B-P®** Bard-Parker®  
SAFETYLOCK™  
Surgical Blade System



**VACUTAINER®**  
SafetyGlide™ Blood  
Collection Assembly




# Your safety. Our concern.

Becton Dickinson, an innovator of quality medical products for more than a century, is also the leader in sharps safety. Our products, educational programs and exposure surveillance systems help to provide you with a safer work environment, whether you are giving injections, drawing blood, or disposing of sharps.

**BD. One Company, Many Solutions.**




**B-D®** Insyte®  
AutoGuard™  
Shielded IV Catheter




**B-D®** SafetyGlide™  
Shielding  
Injection Needle




**B-D®** GUARDIAN®  
Sharps Collectors



**B-D®** Saf-T-Intima™  
IV Catheter  
Safety System



**B-D®** Cannula  
InterLink®  
IV Access System



**VACUTAINER®**  
Safety-Lok™  
Blood Collection Set

**BECTON  
DICKINSON**

For more information, call 1-888-237-2762, press Prompt #1 or visit our web-site at <http://www.bd.com>

InterLink is a trademark of Baxter International Inc.

All other noted trademarks are property of Becton Dickinson and Co., Franklin Lakes, NJ

© 1999 Becton Dickinson and Co.

<https://doi.org/10.1017/S0195941700070156> Published online by Cambridge University Press

# AT LAST

*A hydrogen peroxide  
solution for flexible  
endoscope reprocessing*

*Completely immersed in*

# BENEFITS!

- Premixed formula for automatic and manual reprocessing<sup>1</sup>
- Fast-acting sporicide, tuberculocide, fungicide, virucide, and bactericide<sup>1</sup>
  - High-level disinfection in 30 minutes at room temperature
  - Sterilization in 6 hours at room temperature
- Does not contain glutaraldehyde
  - no aldehyde odor
- Cost-effective reusable solution<sup>1, 2</sup>
- Compatible with metals, plastics, and elastomers such as those in flexible endoscopes<sup>3, 4, 5, 1</sup>

For more information, call  
**1-888-4SPOROX.**

**RECKITT & COLMAN**

The name to trust for your  
infection control needs



# SPOROX™

STERILIZING & DISINFECTING  
SOLUTION

**The simple solution  
for reprocessing**

When your standard of care is high level disinfection. See package insert for complete listings and contraindications.

\*Reusable when used in conjunction with SPOROX test strips. <sup>2</sup>Consult with your endoscope manufacturer before using SPOROX Solution.

References: 1. SPOROX Solution package insert. 2. Data on file, Reckitt & Colman Inc.

Published online by Cambridge University Press

©1997 Reckitt & Colman Inc.