

# 10 good reasons to warm up to a Cold Sterilant

- ① Tuberculocidal in 10 minutes using the EPA quantitative method.
- ② Proven effective against HIV (AIDS virus).
- ③ Pleasant lemon scent.
- ④ Gentle to use, lowest primary dermal irritation score (.34).
- ⑤ Two-year shelf life.
- ⑥ No activator required.
- ⑦ No need to heat to 77°F, effective at room temperature.
- ⑧ Approved 1:4 dilution for an economical .5% high level disinfection.
- ⑨ Won't fog or corrode fiberscopes-safe to use. Contains a corrosion inhibitor.
- ⑩ Longest re-use life of all disinfectants.




The fastest tuberculocidal and longest-lasting sterilizing/disinfecting solution available today.

**WAVICIDE-01**  
(2% GLUTARALDEHYDE)

Helping you turn  
the tide in  
infection control.

➔ Wave Energy Systems, Inc.  
E.P.A. Establishment No. 15136-NJ1  
218 Little Falls Road  
Cedar Grove, N.J. 07009  
800-252-1 125



YOUR  
JOB  
CAN BE  
LESS  
OF A  
THREAT

“HANDS-OFF”


INJECTOR  
**Tubex**<sup>®</sup>

(closed injection system) Ayerst

TUBEX<sup>®</sup>, the unique, *safer*, hands-off injector that puts you farther from the needle and cartridge than any other injection system.

There's no touching...shaking... manipulation of cartridge or needle. Just twist and drop.

Nurses preferred TUBEX almost 3:1 for the prevention of accidental needle stick injuries in a study comparing TUBEX Injector and Carpuject<sup>®</sup>\*1 Which is one more reason to



adopt TUBEX- the injection system designed as if your life depended on it.

**FARTHER FROM  
THE POINT OF  
DANGER**

1 Data on file, Ayerst Laboratories  
\*Carpuject<sup>®</sup> is a registered trademark of Winthrop Pharmaceuticals, Inc

**WYETH-AYERST  
LABORATORIES**

Philadelphia, PA 19101 ©1989, Ayerst Laboratories

5466